<u>küschall</u>® designed for life

OPTIONS & ACCESSORIES

Swiss precision quality.

Welcome to the world of küschall – a world where the user's way of life and aspirations drive everything. Every wheelchair we deliver is unique and built specifically for one individual; always with the aim to exceed the user's expectations.

We strive to ensure our products work perfectly. We work to create solutions designed to make everyday life easy and more convenient. Each wheelchair is hand built and assembled by a specialist fitter in our factory in Switzerland and is then subjected to a stringent quality inspection.

Precision is essential for küschall, and pushes us to bring more and more information to you. With the küschall Options & Accessories brochure, you have all information at hand to see and compare the available parts for our wheelchairs. Thanks to the detailed data it is easier to compare the weight and options to make a safe choice for your customer.

Except any other indication, all weights indicated are based on seat width 400mm and seat depth 400mm with a tolerance value of +/- 15%. Options weights are calculated by pair when requested (e.g. for the wheels, push handles, brakes) and with their related fixations.

All options are subject to compatibility with other options you may have chosen. Always refer to the order form for option availability and compatibility. The reference codes in this brochure are not part numbers, but relate to the order form and are for reference purpose only. Please refer to the electronic spare parts catalogue for part numbers. Merchandising articles items only available as long as on stock.

Contact your nearest küschall specialist for more information.

Follow us on social media











Küschall AG | Benkenstrasse 260 | CH-4108 Witterswil | Switzerland Phone +41 (0) 61 487 70 70 | Fax +41 (0) 61 487 70 79 | www.kuschall.com

UNITED KINGDOM: Invacare Ltd. | Pencoed Technology Park, Pencoed, Bridgend CF35 5AQ Phone +44 (0) 1656 776222 | Fax +44 (0) 1656 776220 | www.invacare.co.uk

IRELAND: Invacare Ltd. | Unit 5, Seatown Business Campus | Seatown Road | Swords, Co. Dublin Phone +353 (0)1 810 7084 | Fax +353 (0)1 810 7085 | www.invacare.ie

Invacare EASTERN EUROPE, MIDDLE EAST, RUSSIA & CIS | Kleiststraße 49 | D-32457 Porta Westfalica Phone +49 (0)57 31 754 540 | Fax +49 (0)57 31 754 541 | edo@invacare.com



Küschall is a registered trademark. Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, Küschall AG reserves the right to modify existing models at any time. Any use of this publication or of parts thereof as well as any reproduction of images must have the written consent of Küschall AG.

2016 Copyright © 09/2016, Küschall AG, Switzerland – All rights reserved

CONTENT

BACKREST UPHOLSTERY	4
MATRX BACKS	5
PUSH HANDLES	6
SEATING	7
FOOTREST	9
SIDEPARTS	10
THE COLOURS	12
FRONT WHEELS	13
REAR WHEELS	14

TYRES	16
AXLE	18
HANDRIMS	19
BRAKES	20
MISCELLANEOUS	22
MERCHANDISING	24
MYKÜSCHALL	28
ADD-ON SOLUTIONS	29



BACKREST UPHOLSTERY











	Product description	Weight	Option Nr.
1	Standard back upholstery, "black premium", tension adjustable	500 gr (including Velcro® straps)	0726
2	Standard back upholstery, black, only for attract family, tension adjustable	523 gr (including Velcro® straps)	0725
3	Aero back upholstery, breathable material, tension adjustable	630 gr (including Velcro® straps)	0735
4	Standard back upholstery, not tension adjustable	223 gr	0700
5	Light back upholstery, not tension adjustable	235 gr	0721

MATRX BACKS









	Product description	Weight *	Option Nr.
1	The Matrx Elite Back provides targeted and comfortable support for enhanced posture and function. Shaped to follow the body's natural contours, the attractive, lightweight design is ideal for all users.	1.2 kg 1.4 kg 1.5 kg	0758 (260 mm) 0759 (310 mm) 0760 (360 mm)
2	The Matrx Elite Deep Back is designed with built-in lateral supports, providing additional trunk support where it is needed most. Hip and scapular cut-outs ensure both clearance and freedom of movement for users seeking a more active lifestyle.	1.6 kg	0756 (310 mm)
3	The Matrx MX1 Back is designed for active individuals who demand peak performance without compromising aesthetic appeal. This latest innovation in back support technology combines superior positioning and function. The full range of hardware options ensure optimal adjustment, while the lightweight carbon fibre shell provides exceptional support and durability.	540 gr 630 gr 770 gr	0753 (230 mm) 0754 (310 mm) 0755 (410 mm)

More information: motionconcepts.com




^{*} Weights are without mounting hardware


PUSH HANDLES













	Product description	Weight*	Option Nr.	
1	Push handles fixed short (73 mm)	332 gr (incl. tubes and fixations)	0751	
2	Push handle fixed long (110 mm)	760 gr (incl. tubes and fixations) 240 gr (excl. tubes and fixations)	0702	
3	Push handles height adjustable integrated	760 gr (incl. tubes and fixations)	0733	
4	Push handles, foldable inactive/active	524 gr (incl. tubes and fixations)	0704	
5	Push handles height adjustable rear seat	740 gr (incl. tubes and fixations)	0703	

^{*} Weight per pair

SEATING















	Product description	Weight	Option Nr.
1	Light seat upholstery, not tension adjustable	143 gr	0615
2	Standard seat upholstery with 1 central velcro, not tension adjustable	200 gr	0610
3	Tension adjustable seat upholstery with velcro surface (soft)	510 gr	0611
4	Seat upholstery bag, small or big, to be fixed under the seat by velcro band (small: 16x34 cm) (large: 24x34 cm)	small: 50 gr large: 65 gr	0616 (small) 0617 (large)
5	Seat bag, click in and removable (dimension: 190x95x40 mm)	100 gr	0606
6	Mobile phone holder, mounted with velcro straps (dimension: 125x70x20 mm)	40 gr	1746

SEATING









	Product description	Weight	Option Nr.
1	Standard seat cushion with visco top (memory foam) or alcantara textile (without visco top)	910 gr	0600
2	The Matrx Vi cushion is designed to provide skin protection and positioning for individuals at moderate to high risk of skin breakdown. The unique three part foam construction targets skin protection where it is needed most, and features excellent stability, support and comfort.	1.4 kg	0663
3	The Matrx Libra cushion uses a combination of high resilient foam and dual layer fluid sac to provide maximum protection around boney prominences and stability to the pelvis. Recommended for all users including those at highest risk of skin breakdown. The Libra Fit Kit components are available to enhance immersion, correct and accommodate orthopedic issues.	2.0 kg	0662



FOOTREST



	Product description	Weight	Option Nr.
1	Flip up footplate, angle adjustable, composite material	280 gr SB38	1500 / 1501
2	Carbon footplate, angle adjustable Carbon flip-up footplate, adjustable, fixed on 1 side when folded	160 gr SB38	1570 / 1515 1566 / 1567
3	Flip up 2-pieces footplate, synthetic material	320 gr SB40	1562
4	2-piece footplate, carbon	240 gr	1578 / 1579
5	Tube footrest, black or raw Aluminum or Titanium	Aluminum: 260 gr Titanium: 210 gr	Aluminum: 1508 Titanium: 1509
6	Footrest cover for tube footrest	100 gr SB42	1513
7	1-piece for rigid, 2-piece for folding wheelchairs	1 pc: 50 gr 2 pc: 100 gr	1543
8	Heel loops	320 gr	1530

SIDEPARTS
















	Product description	Weight*	Option Nr.
1	Carbon clothes guards	360 gr	1400
2	Economy clothes guards, composite material for Compact attract	460 gr	1424
3	Carbon mud guards	320 gr	1405
4	Carbon clothes guards, flip back	200 gr	1401
5	Composite clothes guards, flip back for K-Series attract	320 gr	1424
6	Economy mud guards, composite material removable	600 gr	1425
7	Carbon mud guards, removable, (This option adds 2 cm to the total width of the wheelchair)	800 gr	1408

^{*} Weight per pair




	Product description	Weight*	Option Nr.
1	Siderest height adjustable, carbon, removable, short & long, non-locking	1.2 kg	1416 (short) 1421 (long)
2	Siderest height adjustable, composite, removable, short & long, non-locking	1.6 kg	1416 (short) 1421 (long)
3	Armrest hemi (desk on page 20)	1.1 kg	1422 (right) 1423 (left)
4	Tubular armrest, height adjustable, removable, swing-away, comfort padded	620 gr	1411
5	Armrest, removable, height adjustable, removable polyester cover, rim in Red, Black or Blue (Armpad lenght 280 mm)	980 gr	1445 (black) 1446 (blue) 1447 (red)
6	Pair of armrest covers	214 gr	1455 (black contour) 1456 (blue contour) 1457 (red contour)

6

* Weight per pair

THE COLOURS



FRONT WHEELS
















	Product description	Size	Width	Weight*	Option	Nr.						
1	Starec wheel, available in CPS colours	4" 5"	24 mm 24 mm	_	0809 0819	chrome 0861 0862	techno- anthracite 0803 0813	0808	copper glaze 0830 0833	0807	arctic blue 0831 0832	
2	Super sport	3"	24 mm	180 gr	0800							
3	Comfort and everyday wheel	4" 5"	29 mm 31 mm		Everyday 0801 0810	Comfort 0860 0863						
4	Light wheel	5" 6" 7"	31 mm 32 mm 35 mm	624 gr 788 gr 990 gr	0810 0820 0822							
5	Pneumatic wheel	6″	31 mm	860 gr	0821							
6	FrogLegs wheel Epic Aluminum Wide	4" 5"	36 mm 36 mm	9	black 0796 0798	silver 0797 0799						


REAR WHEELS











15

١		Product description	Weight*	Option Nr.
	1	Küschall Starec, cross spoked Starec hub colourless and rim black - Hub and rim in CPS when selected	22": 1.4 kg 24": 1.5 kg 26": 1.7 kg	0915
	2	High Performance wheels, Starec hub black - or in CPS when selected, black rim	24": 1.4 kg 25": 1.5 kg	0908 (24") 0918 (25")
	3	Spinergy LXK, Hub, rim and spokes black-CPS not applicable if choosen	24": 1.2 kg 25": 1.3 kg	0906 (24") 0912 (25") NEW
	4	Wheels with integrated hub brakes, 24 silver spokes, black hub and rim - CPS not applicable if choosen	24": 1.6 kg	0927
	5	Mountain Bike wheels, complete wheel with Mountain Bike tyre	24": 4.0 kg	0913
	6	Active wheels, complete wheels: light tyre 7,5 bar and alu grey handrims with 6 points of fixation (screws)	24": 3.1 kg	0921
	7	Drum brake wheel, automatically delivered with drum brake	24": 3.5 kg	1304
		One arm drive system, left or right side (with light wheel 7 bar tyres, alu handrims)	24": 5.2 kg	1118 (right) / 1119 (left)
		Anti flat tyre tube protection (placed between tyre and rim)	110 gr (1 pc)	1020
	10	1: High Performance / Spinergy LXK / Active 2: Starec	1: 30 gr 2: 10 gr	1710 (mounted on spokes) 1711 (enclosed)
	11	Emotion adapter plate / Twion adapter plate	1.2 kg	0916 (E-Motion) / 0926 (Twion)

* Weight per pair


TYRES



	Product description	Weight*	Option Nr.
1	Mountain bike tyre, only available as complete wheel (page 13)	24": 750 gr (without wheel)	0913
2	Marathon Plus tyre with Smartguard, 5mm layer of protective rubber infill, black	24": 1.1 kg	1010
3	Profile tyre 7 bar, grey	24": 600 gr	1002
4	Solid tyre grey	24": 1.0 kg	1007
5	Solid KIK tyre, black	24": 900 gr	1006
6	Rear wheel solid tyre grey	24": 1.1 kg	1004
7	Schwalbe one tyre, V guard puncture resistant	24": 500 gr	1011
8	Schwalbe rightrun tyre, black/blue	24": 600 gr	1013
9	Schwalbe rightrun tyre, black/grey	24": 600 gr	1014



AXLE



	Product description	Weight	Option Nr.
1	Aluminum axle black	320 gr	0277
2	Carbon axle	220 gr	0278
3	Vario-ax, adjustable camber 1°,3°,7°,10°	600 gr	0276
4	Rear wheel extension, 3 or 5 positions	198 gr / 270 gr	0248
5	Axle protection for carbon axles	62 gr	1750
6	Rear wheel extension, for folding wheelchairs	296 gr	0249
7	Quick release axle, available in CPS colours	100 gr	N/A
8	Tetra clip	62 gr	0909


HANDRIMS




^{*} Weight per pair


BRAKES


















	Product description	Weight*	Option Nr.
1	Active brake	290 gr	1303
	Performance brake	420 gr	1302
	Push brakes lever straight	380 gr	1301
	Push brakes, bent	360 gr	1300
5	Pull brake	400 gr	1306
6	Brake lever extension	N/A	1310
	Push brakes, lever to fold down, only for "attract" family	10 gr	1312

^{*} Weight per pair

MISCELLANEOUS





















		Product description	Weight	Option Nr.
-		Active antitip	510 gr 1pc	1705
4	2	Anti tipper	520 gr 1pc	1702 (right), 1703 (left)
3	3	Tipper aid	114 gr 1pc	1708 (right), 1709 (left), 1718 (both sides)
4		Cane holder	154 gr	1706 (right), 1707 (left), 1719 (both sides)
į	5	Transit wheels	1.1 kg the pair	1700
•	6	Rear frame protection	20 gr the pair	1729
-	7	Air man	1.9 kg	1716
8	8	Air pump	66 gr	1715
9	9	Bottle holder with küschall bottle	140 gr (including the bottle)	1745
1	.0	Desk, siderest on page 9 required	1.5 kg	1760

















	Product description	Weight	Option Nr.
1	Luggage holder, complete pair, active/inactive	400 gr the pair	1744
2	Front frame protectors (2 pieces)	80 gr the pair	1739
3	Lateral support	512 gr	N/A
4	Travel bag for küschall champion	Small: 2.0 kg Big: 2.8 kg	1740
5	Wheelbag, for rear wheels, fits all size wheels	778 gr	1728
6	Transportation kit including posture belt	350 gr	1770
7	Rear hight adjustment set for küschall K-Series	40 gr the pair	1727
8	Suspension for küschall K-Series, max. 100kg user weight	770 gr	0630 (40-65kg) 0631 (65-85kg) 0632 (85-100kg)

MERCHANDISING


















	Product description	Option Nr.
1	T-shirt, cotton, black, "crystals swooshes"	women: S: 1517777, M: 1517778, L: 1517779, XL: 1517780
2	Scarf, silk, black, "crystal swooshes"	unisex: 1517781
3	T-shirt, cotton, black, "küschall Team"	women: S: 1531189, M: 1531379, L: 1531380, XL: 1531381 men: M: 1531382, L: 1531383, XL: 1531463, XXL: 1531464
4	Sweat-Shirt, cotton, black, "küschall"	unisex: S: 1550289, M: 1550290, L: 1550291, XL: 1550292
5	T-shirt, cotton, black/yellow, "küschall Champion"	women: M: 1517773,L: 1517774 men: S: 1517728, M: 1517729, XL: 1517731, XXL: 1517732
6	T-Shirt, cotton, black/red, "R33"	unisex: S: 1498701, M: 1498702, L: 1498734, XL: 1498735, XXL: 1498776















MERCHANDISING



















-		Product description	Option Nr.
	1	Carbon watch, "küschall"	unisex: 1498740
	2	USB-stick, 2GB, "küschall"	1513362
	3	Pen-set, ballpoint & pencil, "küschall"	1486577
	4	3-Pen-set, red, white, black, "küschall"	1486578
	5	Colour samples, Aluminum, frame & CPS. Availability of the colours "Caramel Bronze, Le Corbusier red, Neon: green,red,yellow,pink" on request	1531815 Standard colours & CPS without the 4 attract colours: 1531815, 4 attract colours: 1584879
	6	POS box, cardboard, including frame and handrim samples, without catalogues	1561019
	7	Wheelchair backpack	unisex: 1532047
	8	Paper bag R33	on request
	9	Küschall club display, without cards	1605587















	Product description	Option Nr.
1	MYküschall clips, set of 50 pcs	1578806
2	MYküschall key ring, 1 pc	1578805
3	MYküschall pencil, set of 4 pcs	1578807
4	Measuring tape, milimetres, "küschall"	1486620
5	Weight scale, metal, "The KSL"	1525328
6	Spirit level, acrylic resin, "küschall"	1567762
7	Tool Kit, "küschall"	1717

MYKÜSCHALL

YOU ARE AN ESSENTIAL PART OF THE MYKÜSCHALL TEAM: YOU HAVE THE IDEA, WE TURN IT INTO REALITY!

Bringing together more than 35 years of innovation and experience, MYküschall is a dedicated service, which creates and offers «designed for life» customized wheelchairs. From technical to aesthetic design customizations, we can optimize your wheelchair to your needs within a wide range of possibilities.

Based on the therapeutic guidelines that you give us, we can customize your next wheelchair so that it meets your postural requirements, thereby improving your comfort and the overall driving performance. From frame colours to special upholstery for seat and backrests, there are plenty of customizations possible.















- Customized front frameDifferent tyres
- 3 Coloured Spinergy LXK spokes
- 4 Individual handrim colour
- 5 Individual frame colour
- 6 Coloured Spinergy LXK spokes

Your satisfaction is key for us! We have brought togethe experienced küschall employees who have all worked at the Swiss factory for many years. From customer service to engineering, manufacturing and quality, these dedicated team members work together to provide safe solutions to your customization requests.

ADD-ON SOLUTIONS

With its wheelchair add-on drives Alber offers a variety of products allowing to expand the function of a manual wheelchair for a better adaption to the needs of the wheelchair users. All Alber products can be mounted onto your Küschall wheelchair with a special developed bracket kit which will not change the wheelchair characteristics nor its functionality.



TWION

Efficient, smart and fun to drive, that's what the lightest and fastest power add-on twion stands for. The silent direct drive ensures high agility and simple handling allowing to speed up to 10km/h. The wheels only weigh 6kg including batteries. With the optionally available Mobility Plus package, the user can use even more features in the twion App and communicate directly with the twion wheels via the smartphone.



E-MOTION

The Alber e-motion merges active mobility with therapeutic benefit. Even wheelchair users with impaired strength can move around independently with the help of e-motion. Using the remote control the e-motion can be switched on and off from an upright sitting position and in addition, two selectable power stages assure optimal driving convenience – precisely adapted to each user's personal mobility restriction and needs.



E-FIX

Small, light and swift. The e-fix converts a manual wheelchair into a lightweight electric wheelchair. It's so easy to control your wheelchair with e-fix you can get used to it almost at once offering excellent comfort thanks to the compact control unit. The clear, backlit colour display shows all the important information at a single glance and the two quiet but powerful drive units allow turning on the spot and manoeuvring even in confined spaces, or safely and powerfully negotiating slopes. The small Li-lon battery pack can be handled easily with only 2kg in weight and a width of 8cm



More information: www alber de

NOTES



